

Exhibit V

Internal Method Comparison Between EPA 300 and SM4500 Cl B

Sample	Chloride mg/kg via EPA300	Chloride mg/kg via SM4500	RPD	BS/BSD REC 300	BS/BSD REC 4500
1	<10	<16	0.00	92.7/93.6	104/104
2	291	256	12.80	92.7/93.6	104/104
3	<50	32	0.00	92.7/93.6	104/104
4	648	576	11.76	92.7/93.6	104/104
5	1190	1376	-14.50	92.7/93.6	104/104
6	13.4	16	-17.69	93.6/94.6	104/108
7	15.2	<16	0.00	93.6/94.6	104/108
8	27300	26800	1.85	93.6/94.6	104/108
9	1080	1088	-0.74	95.5/95.6	108/108
10	574	624	-8.35	95.5/95.6	108/108
11	701	656	6.63	95.5/95.6	108/108
12	53.3	64	-18.24	95.5/95.6	108/108
13	6250	6000	4.08	95.5/95.6	108/108
14	67.5	64	5.32	94.3/93.8	104/108
15	68.3	64	6.50	94.3/93.8	104/108
16	243	272	-11.26	94.3/93.8	104/108
17	108	128	-16.95	94.3/93.8	104/108
18	92.2	96	-4.04	94.3/93.8	104/108
19	103	128	-21.65	94.3/93.8	108/104
20	227	240	-5.57	92.8/93.0	100/108
21	306	352	-13.98	92.8/93.0	100/108
22	455	544	-17.82	92.8/93.0	100/108
23	500	528	-5.45	92.8/93.0	100/108
24	218	272	-22.04	92.8/93.0	100/108
25	2760	2599	6.01	97.7/98.5	100/108
26	474	512	-7.71	97.7/98.5	100/108
27	495	544	-9.43	92.3/91.8	112/108
28	920	816	11.98	92.3/91.8	112/104
29	292	320	-9.15	92.3/91.8	104/104
30	211	256	-19.27	93/93.3	104/108
31	4030	4320	-6.95	90.9/91.4	108/104
32	2120	2520	-17.24	90.9/91.4	108/104
33	1610	1630	-1.23	90.9/91.4	108/104
34	4740	4560	3.87	90.9/91.4	108/104
35	850	944	-10.48	90.9/91.4	108/104
36	758	848	-11.21	90.9/91.4	108/104

%RPD = Relative Percent Difference

BS = Blank Spike

BSD = Blank Spike Duplicate

%REC = Percent Recovery

EPA 300 analysis by Green Analytical Laboratories. SM4500 Cl B analysis by Cardinal Laboratories.