

ENVIROTECH LABS

PRACTICAL SOLUTIONS FOR A BETTER TOMORROW

EPA METHOD 8040 PHENOLS Quality Assurance Report

Client:	QA/QC	Project #:	N/A
Sample ID:	Matrix Duplicate	Date Reported:	07-25-00
Laboratory Number:	H805	Date Sampled:	N/A
Sample Matrix:	Water	Date Received:	N/A
Preservative:	Cool	Date Extracted:	N/A
Condition:	Cool & Intact	Date Analyzed:	07-24-00
		Analysis Requested:	TCLP

Parameter	Sample Result (mg/L)	Duplicate Result (mg/L)	Detection Limit (mg/L)	Percent Difference
o-Cresol	ND	ND	0.020	0.0%
p,m-Cresol	ND	ND	0.040	0.0%
2,4,6-Trichlorophenol	ND	ND	0.020	0.0%
2,4,5-Trichlorophenol	ND	ND	0.020	0.0%
Pentachlorophenol	ND	ND	0.020	0.0%

ND - Parameter not detected at the stated detection limit.

QA/QC Acceptance Criteria:	Parameter	Maximum Difference
	8040 Compounds	30.0%


References: Method 1311, Toxicity Characteristic Leaching Procedure Test Methods for Evaluating Solid Waste, SW-846, USEPA, July 1992.


Method 3510, Separatory Funnel Liquid-Liquid Extraction, Test Methods for Evaluating Solid Waste, SW-846, USEPA, July 1992.

Method 8040, Phenols, Test Methods for Evaluating Solid Waste, SW-846, USEPA, Sept. 1986.

Note: Regulatory Limits based on 40 CFR part 261 subpart C section 261.24, July 1, 1992.

Comments: QA/QC for sample H805.


Analyst


Review