

Gauge No. 1496		Depth 12619'		Clock No. 10567		24 hour		Ticket No. 604312	
First Flow Period		First Closed In Pressure		Second Flow Period		Second Closed In Pressure		Third Flow Period	
Time Defl. .000"	PSIG Temp. Corr.	Time Defl. .000"	PSIG Temp. Corr.	Time Defl. .000"	PSIG Temp. Corr.	Time Defl. .000"	PSIG Temp. Corr.	Time Defl. .000"	PSIG Temp. Corr.
0 .0000	4167.4	.0000	1033.1	.0000	997.7	.0000	1121.6		
1 .0867*	3316.2	.0136**	1159.2	.0536***	1028.7	.0338***	1493.3		
2 .1170-0	3042.3	.0443	1519.9	.1038	1055.3	.0744	2026.6		
3 .1534	2139.6	.0749	1913.7	.1541	1073.0	.1150	2580.9		
4 .2201	1413.7	.1055	2354.7	.2043	1090.7	.1556	2960.0		
5 .2868	1115.0	.1362	2745.0	.2545	1103.9	.1961	3198.2		
6 .3535	1030.9	.1668	3026.7	.3048	1112.8	.2367	3347.4		
7 .4200	1033.1	.1975	3231.6	.3550	1121.6	.2773	3452.1		
8		.2281	3369.7			.3179	3530.0		
9		.2587	3469.9			.3585	3590.2		
10		.2894	3545.6			.3990	3630.2		
11		.3200	3601.3			.4396	3665.9		
12						.4802	3692.6		
13						.5208	3717.1		
14						.5614	3737.2		
15						.6020	3755.0		

Gauge No. 1516		Depth 12922'		Clock No. 5991		24 hour	
First Flow Period		First Closed In Pressure		Second Flow Period		Second Closed In Pressure	
Time Defl. .000"	PSIG Temp. Corr.	Time Defl. .000"	PSIG Temp. Corr.	Time Defl. .000"	PSIG Temp. Corr.	Time Defl. .000"	PSIG Temp. Corr.
0 .0000	4263.6	.0000	1148.7	.0000	1107.7	.0000	1237.0
1 .0852*	3421.7	.0131**	1297.4	.0515***	1153.0	.0329***	1637.9
2 .1150-0	3145.6	.0426	1670.2	.1018	1176.7	.0723	2194.8
3 .1508	2242.4	.0721	2090.9	.1510	1193.9	.1117	2718.6
4 .2163	1517.2	.1016	2510.8	.2003	1206.8	.1512	3091.3
5 .2819	1226.2	.1311	2891.7	.2495	1219.8	.1906	3321.7
6 .3475	1142.2	.1606	3167.3	.2988	1232.7	.2301	3471.7
7 .4130	1148.7	.1900	3360.8	.3480	1237.0	.2695	3573.9
8		.2195	3497.8			.3089	3647.8
9		.2490	3593.4			.3484	3704.3
10		.2785	3665.2			.3878	3747.8
11		.3080	3717.3			.4273	3782.6
12						.4667	3810.8
13						.5061	3832.6
14						.5456	3854.3
15						.5850	3871.7

Reading Interval 20

REMARKS: *First interval is equal to 26 minutes. ** = 4 minutes. *** = 16 minutes. **** = 10 minutes C = choke change