

Continued

Page 2 - Lorene 1-1 C-103 - tubing & perforations

- 01-20-98 Load & test w/44 bbls 2% KCL, treat perfs 2' spacing on PPI tool @11,619'-11,576', total 245 BLWTR, avg rate .5 BPM, max rate .5BPM, avg PSI 500#, max PSI 2500#, ISIP 0#, release pkr, RU slickline, TIH, fish ICV, SDFN.
- 01-21-98 TOOH w/tbg & PPI tool, PU & TIH w/5-1/2" model R Pkr, 2-1/4" SN & 372 jts 2-7/8" tbg, set pkr @10,510', ND BOPs, NU wellhead,