

CAVINESS FEDERAL NO. 1
CONT.

KC1 WATER CONTAINING 2 GPT OF TFA-380B OR EQUIVALENT.
POH WITH TUBING. RDMO STIMULATION COMPANY

RIH WITH SN, TUBING ANCHOR, AND 2 7/8" TUBING. SET ANCHOR AT PREVIOUSLY SET DEPTH
ND BOP. NU WH. RIH WITH PMP AND RODS. SPACE OUT PMP AND CONNECT TO PUMPING UNIT.
RDMO PULLING UNIT.