

DST SUMMARY

8-21-80 Test: Bone Spring (9598'-9626') - 17 min. preflow; open good blow gas to surface. ISI 1 hr. Final flow tool open 1 hr., oil to surface in 23 min. flowing press. @ end of period, 450# on 1/2" surface CK. Est. 50 BBL oil to pit. FSI 2 hr. IH 4658, 15 min. preflow 1059-1799, ISI 3852, FF 1 hr. 1820-4434. FSI 3852, FH 4700. Reverse out 75 BBL heavy gas cut oil to pits. BH Sampler 1400# contains 1100 cc oil, no wtr. 2 CF gas. BHT 148°.