

Form C-103

L. Van Etten Well No. 2

Page 2

- 15) Set RTTS packer at 3642' and acidized perfs. 3676-96' with 1250 gallons 15% NE Acid. Flushed with 14 bbls. water. Reset packer at 3680' and tested perfs. 3685-96'.
- 16) Reset packer at 3642' and treated perfs. 3676-96' with 3000 gallons 15% RA-1 acid and 6 ball sealers. Flushed with 35 bbls. water.
- 17) Ran 119 jts. (3659') of 2-7/8" OD tubing and set packer at 3659', seating nipple 3653-3654'.
- 18) Tested.