

South Eunice Unit No. 43

Pg 3

11. Load backside w/ mixture of 7 gals Conoco No. 1 pkr fluid plus 4 gals Conoco No. 2 pkr fluid in 70 bbls 2% KCH TFW. Hatch Abq into pkr at 3550'.
12. Acidize interval 3680'-3800' at 4-6 BPM w/ max surface pressure of 2500 psi using 120 bbls 15% HCl-NE-FE acid & 400# graded rock salt mixed in 3 bbls gelled 10# Irone in 3 stages. Flush w/ 100 bbls 2% KCH TFW w/ 1 gal Adomall / 1000 gals. Swab well.
13. Return well to injection. Let well stabilize at least 1 week.
14. Run injection profile survey.