

BRUNSON ARGO #13, 6580 TD.

Blinebry & Tubb: P 0 BO + 8 BLW/24 hrs, 82 MCFGPD. - 1466 BLW.

BRUNSON-ARGO #13, 6580 TD.

Blinebry & Tubb: P 0 BO + 8 BLW/24 hrs, 83 MCFGPD. - 1458 BLW.

BRUNSON-ARGO #13, 6580 TD.

Blinebry & Tubb: P 0 BO + 8 BLW/24 hrs, 88 MCFGPD. - 1450 BLW.

BRUNSON ARGO #13, 6580 TD.

Blinebry & Tubb: P 0 BO + 4 BLW/24 hrs, 92 MCFGPD. -1446 BLW.

BRUNSON ARGO #13, 6580 TD.

Blinebry & Tubb: P 0 BO + 24 BLW/3 days, 95 MCFGPD. - 1422 BLW.

BRUNSON ARGO #13, 6580 TD.

Blinebry & Tubb: P 0 BO + 8 BLW/24 hrs, 100 MCFGPD. - 1414 BLW.

BRUNSON ARGO #13, 6580 TD

Blinebry & Tubb: P 0 BNO + 10 BLW/ 24 hrs, 74 MCF/D. -1404 BLW.

BRUNSON ARGO #13, 6580 TD.

Blinebry & Tubb: P 0 BO + 8 BLW/24 hrs, 74 MCFGPD. - 1396 BLW.

BRUNSON ARGO #13, 6580 TD.

Blinebry & Tubb: P 0 BO + 18 BLW/24 hrs, 74 MCFGPD. - 1378 BLW.

BRUNSON ARGO #13, 6580 TD.

Blinebry & Tubb: P 0 BO + 27 BLW/3 days, 78 MCFGPD.
- 1351 BLW.

BRUNSON ARGO #13, 6580 TD.

Blinebry & Tubb: P 0 BO + 8 BLW/24 hrs, 78 MCFGPD. - 1343 BLW.

BRUNSON ARGO #13, 6580 TD.

Blinebry & Tubb: P 0 BO + 8 BLW/24 hrs, 74 MCFGPD. - 1335 BLW.

BRUNSON ARGO #13, 6580 TD.

Blinebry & Tubb: P 0 BO + 10 BLW/24 hrs, 65 MCFGPD. - 1325 BLW.

BRUNSON ARGO #13, 6580 TD.

Blinebry: Eunice Well Serv MIRU PU, pulled rods and pump,
set thg plug in "F" nipple @ 5899, prep to POP
and test Blinebry. - 594 BLW Est.

(Tubb: Shut off. - 731 BLW Est.)

BRUNSON ARGO #13, 6580 TD.

Blinebry: Ran rods and pump,
P 0 BO + 43 BLW + 220 MCFGPD/3 days. - 551 BLW.

Tubb: SI. - 731 BLW.

BRUNSON ARGO #13, 6580 TD.

Blinebry: P 0 BNO + 28 BLW + 215 MCFG/24 hrs. - 523 BLW.

Tubb: SI. - 731 BLW.

BRUNSON ARGO #13, 6580 TD.

Blinebry: P 0 BO + 11 BLW + 225 MCFG/24 hrs. - 512 BLW.

Tubb: SI. - 731 BLW.